

BMT CTN 1202
Research Sample Types and Collection Schedule

Sample Type	Sample Aliquot Volume & Quantity	Time Point							
		Day -1/0	Day 7	Day 14	Day 21	Day 28	Day 42	Day 56	Day 90
Serum	0.25 mL Max. 10 aliquots	X	X	X	X	X	X	X	X
Plasma (EDTA)	0.25 mL Max. 10 aliquots	X	X	X	X	X	X	X	X
PAXgene Lysates	2.5 mL-fill frozen PAXgene tubes Max. 2 aliquots				X			X	X
Buffy Coat Cells	Buffy coat cells isolated from 5 mL tube submitted for plasma processing. Single aliquot of cells for future genomic DNA isolation		X	X	X	X	X	X	X

Note: Peripheral blood samples are shipped overnight to the BMT CTN Central Processing Laboratory/Repository for processing and sample aliquot storage. Plasma, serum, buffy coat cells and PAXgene sample aliquots are being stored at -80°C.