

**BMT CTN 0802**  
**Research Sample Types and Collection Schedule**

Sample Type	Sample Aliquot Volume & Quantity	Pre-Treatment Baseline	Post-treatment	
			Day 28	Day 56
Whole Blood (source donor DNA)	1 mL 3 aliquots	X		
Buccal Cells (source patient DNA)	3-4 buccal swabs	X		
Plasma (EDTA)	0.5 mL 2-12 aliquots	X	X	X

**Note:** Sample aliquot processing being performed by the clinical sites. Frozen sample aliquots are periodically batch-shipped to the BMT CTN Repository. All samples being stored at -80°C.